

REMOVE A WORKER

VAULT



[DISCLAIMER]

This presentation document has been prepared by Vault Intelligence Limited ("Vault") and is intended for off line demonstration, presentation and educational purposes.

The information is not version controlled and as such is subject to change without notice.

Vault does not accept any responsibility or obligation to inform users of this presentation of such changes.

This includes any copies of this presentation taken and modified outside of the Vault domain.

[Remove a Worker] Removal Methods



There are two ways to **remove** a Worker from the *Worker Register*:

Archive Worker. This removes the Worker from the Worker Register, but keeps a copy of their data should it be required in the future. This is the recommended method of removing a Worker

Delete the Worker. This removes the Worker from the *Worker Register* and deletes all data relating to the Worker. **Note**: Your organisation may have disabled this feature

In this guide we will focus the **Archive** method


ARCHIVE WORKER

[Remove a Worker] Worker Register



Archiving a Worker is done via the *Worker Register*, this is found under *People* in the *menu* bar


The Worker Register shows all current Workers within the organisation, along with key information about each individual Worker



[Remove a Worker] Final Medical



If a Final Medical Examination has been carried out for the Worker leaving, **update** this before archiving the Worker record




[Remove a Worker] Archive Worker



Locate the Worker in the Worker Register and click the Actions button, then select Archive. If you're sure, click Yes to confirm the action, when prompted

The Worker record is now archived, meaning the Worker has been removed from the *Worker Register*, but data has been safely stored in Vault



[Remove a Worker] Unarchive Worker




If you've made a mistake or you need to **restore** a Worker for any reason, it's easy to **unarchive** a Worker.

To locate an archived Worker:

Click Show Archived People in the Filters sidebar

Supply some of the Workers detail as a search term (e.g. Name, Employee No.)

Then, **click** Apply



[Remove a Worker] Unarchive Worker cont.



With the Archived Worker now located, **click** *Actions*, then *Unarchive*

The Worker record is now **restored** to the *Worker Register*

